



## Product specification

**lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, suspension cells**



## Product information

Product number	94.6077.333
----------------	-------------

## Product descriptions

Name	lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, suspension cells
------	---

Product number	94.6077.333
----------------	-------------

Product description	lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, for suspension cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag
---------------------	---

## Material & colours

Frame material	Polystyrene (PS)
----------------	------------------

## Size

Diameter	35 mm
----------	-------

Height of product	6 mm
-------------------	------

Foil thickness	25 µm
----------------	-------

## Product characteristics

Growth area	6.35 cm <sup>2</sup>
-------------	----------------------

Design	with foil base
--------	----------------

Base shape	Flat base
Vessel type	Dish
Closure type	Slip-on lid
Cell type	suspension cells

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	15000
Height of case	265 mm
Width of case	253 mm
Depth of case	348 mm
Volume of case	0.0233 cbm
Weight of case	1.91 kg
EAN inner box	4038917255141
EAN case	4038917255134

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT