



## Product specification

**lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, suspension cells**



## Product information

|                |             |
|----------------|-------------|
| Product number | 94.6077.333 |
|----------------|-------------|

## Product descriptions

|      |   |
|------|---|
| Name | lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, suspension cells |
|------|---|

|                |             |
|----------------|-------------|
| Product number | 94.6077.333 |
|----------------|-------------|

|                     |   |
|---------------------|---|
| Product description | lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, for suspension cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag |
|---------------------|---|

## Material & colours

|                  |       |
|------------------|-------|
| Product material | lumox |
|------------------|-------|

|                |                  |
|----------------|------------------|
| Frame material | Polystyrene (PS) |
|----------------|------------------|

## Size

|          |       |
|----------|-------|
| Diameter | 35 mm |
|----------|-------|

|                   |      |
|-------------------|------|
| Height of product | 6 mm |
|-------------------|------|

|                |       |
|----------------|-------|
| Foil thickness | 25 µm |
|----------------|-------|

## Product characteristics

|             |                      |
|-------------|----------------------|
| Growth area | 6.35 cm <sup>2</sup> |
|-------------|----------------------|

|              |                  |
|--------------|------------------|
| Design       | with foil base   |
| Base shape   | Flat base        |
| Vessel type  | Dish             |
| Closure type | Slip-on lid      |
| Cell type    | suspension cells |

## Purity & certification

|                  |  |
|------------------|--|
| Product category | no medical device   IVD                              |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Electron irradiation                                 |
| Pyrogen-free     | yes  |
| Batched          | yes  |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 50            |
| Piece(s) / inner box   | 50            |
| Piece(s) / case        | 250           |
| Piece(s) / pallet      | 15000         |
| Height of case         | 265 mm        |
| Width of case          | 253 mm        |
| Depth of case          | 348 mm        |
| Volume of case         | 0.0233 cbm    |
| Weight of case         | 1.91 kg       |
| EAN inner box          | 4038917255141 |
| EAN case               | 4038917255134 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT