



Product specification

Centrifuge, SMC 6 V, supply voltage: 115 V, swing-out rotor 6-fold



Product information

Product number	90.184.401
----------------	------------

Product descriptions

Name	Centrifuge, SMC 6 V, supply voltage: 115 V, swing-out rotor 6-fold
Product number	90.184.401
Product description	Centrifuge, SMC 6 V, supply voltage: 115 V, incl. swing-out rotor 6-fold and tube holder 1-fold for tubes up to 100 mm length, for separation of serum/plasma from the blood cells, swing-out rotor 6-fold

Product characteristics

Design	swing-out rotor 6-fold
Supply voltage	115 V

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / pallet	36
Weight of case	5.8 kg
EAN inner box	4038917534420
EAN case	4038917534413

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT