



Product specification

Automatic-Sarpette® US, Pipetting aids



Product information

Product number	90.189.203
----------------	------------

Product descriptions

Name	Automatic-Sarpette® US, Pipetting aids
Product number	90.189.203
Product description	Automatic-Sarpette® US, Pipetting aids, scope of delivery: charging cable, table stand, wall holder, replacement filters 0.45 µm and 0.2 µm, white/red, ABS, 1 piece(s)/case

Material & colours

Product material	Polypropylene (PP)
Color of product	white/red

Size

Width of product	135 mm
Length of product	40 mm
Height of product	170 mm

Product characteristics

Battery type	Nickel metal hydride battery
--------------	------------------------------

Purity & certification

Product category	Appliance, according to machinery or equipment directive, CE
CE certified	CE - manufacturer's self-declaration
Batched	no

Packaging

Piece(s) / pallet	112
Height of inner box	125 mm
Width of inner box	303 mm
Depth of inner box	333 mm
Height of case	125 mm
Width of case	303 mm
Depth of case	333 mm
Volume of case	0.0126 cbm
Weight of case	1.36 kg
EAN inner box	4038917418966
EAN case	4038917418959

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT