



## Product specification

**Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for**



## Product information

Product number	92.180.611
----------------	------------

## Product descriptions

Name	Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M 2000
Product number	92.180.611
Product description	Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M 2000, 1 piece(s)/case

## Product characteristics

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / case	10
Piece(s) / pallet	10

Volume of case	0.0122 cbm
Weight of case	2.21 kg
EAN inner box	4038917101998
EAN case	4038917075756

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT