



## Product specification

**lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, Adherent cells**



## Product information

Product number	94.6077.410
----------------	-------------

## Product descriptions

Name	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, adherent cells
Product number	94.6077.410
Product description	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, for adherent cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

## Material & colours

Product material	lumox
Frame material	Polystyrene (PS)

## Size

Diameter	50 mm
Height of product	12 mm
Foil thickness	25 µm

## Product characteristics

Growth area	20.35 cm²
-------------	-----------

Design	with foil base
Base shape	Flat base
Vessel type	Dish
Closure type	Slip-on lid
Cell type	adherent cells

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	50
Piece(s) / inner box	50
Piece(s) / case	200
Piece(s) / pallet	7200
Height of case	125 mm
Width of case	162 mm
Depth of case	362 mm
Volume of case	0.0073 cbm
Weight of case	4.59 kg
EAN inner box	4038917255233
EAN case	4038917255226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT