



Product specification

**Block rotor, for 40 tubes up to 11.5 mm diameter, for
SARMIX® M 2000**



Product information

Product number	92.180.612
----------------	------------

Product descriptions

Name	Block rotor, for 40 tubes up to 11.5 mm diameter, for SARMIX® M 2000
Product number	92.180.612
Product description	Block rotor, for 40 tubes up to 11.5 mm diameter, for SARMIX® M 2000, 1 piece(s)/case

Product characteristics

Storage capacity	40
------------------	----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	10
Piece(s) / pallet	10
Height of case	150 mm

Width of case	240 mm
Depth of case	350 mm
Volume of case	0.0126 cbm
Weight of case	1.71 kg
EAN inner box	4038917102001
EAN case	4038917075763

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT