



## Product specification

**Block rotor, for 40 tubes up to 8.5 mm Ø (S-Sedivette®),  
for SARMIX® M 2000**



## Product information

Product number	92.180.615
----------------	------------

## Product descriptions

Name	Block rotor, for 40 tubes up to 8.5 mm Ø (S-Sedivette®), for SARMIX® M 2000
Product number	92.180.615
Product description	Block rotor, for 40 tubes up to 8.5 mm Ø (S-Sedivette®), for SARMIX® M 2000, 1 piece(s)/case

## Product characteristics

Storage capacity	40
------------------	----

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / case	10
Piece(s) / pallet	10
Volume of case	0.0122 cbm

Weight of case	2.71 kg
EAN inner box	4038917102032
EAN case	4038917074315

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT