



## Product specification

**Block rotor, for 14 tubes up to 28 mm diameter, for  
SARMIX® M 2000**



## Product information

Product number	92.180.616
----------------	------------

## Product descriptions

Name	Block rotor, for 14 tubes up to 28 mm diameter, for SARMIX® M 2000
Product number	92.180.616
Product description	Block rotor, for 14 tubes up to 28 mm diameter, for SARMIX® M 2000, 1 piece(s)/case

## Product characteristics

Storage capacity	14
------------------	----

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / case	10
Piece(s) / pallet	10
Height of case	990 mm

Width of case	990 mm
Depth of case	990 mm
Volume of case	0.0122 cbm
Weight of case	3.01 kg
EAN inner box	4038917160858
EAN case	4038917160841

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT