



Product specification

Block rotor, for 24 tubes up to 17 mm Ø (15 ml tubes), for SARMIX® M 2000



Product information

Product number	92.180.617
----------------	------------

Product descriptions

Name	Block rotor, for 24 tubes up to 17 mm Ø (15 ml tubes), for SARMIX® M 2000
Product number	92.180.617
Product description	Block rotor, for 24 tubes up to 17 mm Ø (15 ml tubes), for SARMIX® M 2000, 1 piece(s)/case

Product characteristics

Storage capacity	24
------------------	----

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	10
Piece(s) / pallet	10
Volume of case	0.0122 cbm
Weight of case	2.31 kg

EAN inner box	4038917245906
EAN case	4038917245890

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT