



Product specification

Tube holder, SMC 6, for tubes up to 100 mm



Product information

Product number	92.184.551
----------------	------------

Product descriptions

Name	Tube holder, SMC 6, for tubes up to 100 mm
Product number	92.184.551
Product description	Tube holder, SMC 6, for tubes up to 100 mm

Purity & certification

Pyrogen-free	no
--------------	----

Packaging

Piece(s) / pallet	1944
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.005 cbm
Weight of case	0.22 kg
EAN inner box	4038917162920
EAN case	4038917162890

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT