



Product specification

Cap, LC 6, natural, for Tube holder



Product information

Product number	92.184.555
----------------	------------

Product descriptions

Name	Cap, LC 6, natural, for Tube holder
Product number	92.184.555
Product description	Cap, LC 6, natural, for Tube holder

Material & colours

Color of product	natural
------------------	---------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	60
Piece(s) / pallet	19440
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.43 kg

EAN inner box

4038917102483

EAN case

4038917073257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT