



Product specification

Tube holder, LC 24, red, for 50 ml centrifuge tubes



Product information

Product number	92.184.752
----------------	------------

Product descriptions

Name	Tube holder, LC 24, red, for 50 ml centrifuge tubes
Product number	92.184.752
Product description	Tube holder, LC 24, red, for 50 ml centrifuge tubes

Material & colours

Color of product	red
------------------	-----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / pallet	324
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm

Weight of case	0.54 kg
EAN inner box	4038917102537
EAN case	4038917073288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT