




## Product specification

**lumox® multiwell, Cell culture plate, with foil base, 96 well, 20 piece(s)**



## Product information

Product number	94.6120.096
----------------	-------------

## Product descriptions

Name	lumox® multiwell, Cell culture plate, with foil base, 96 well, 20 piece(s)
------	--

Product number	94.6120.096
----------------	-------------

Product description	lumox® multiwell, Cell culture plate, with foil base, 96 well, frame: black, frame: PS, for adherent cells, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister
---------------------	---

## Material & colours

Frame material	Polystyrene (PS)
Color of frame	black

## Size

Width of product	85.48 mm
Length of product	127.76 mm
Height of product	17.2 mm
Foil thickness	50 µm

## Product characteristics

Growth area	0.34 cm <sup>2</sup>
Number of wells	96
Design	with foil base
Base shape	Flat base
Vessel type	Plate
Closure type	Slip-on lid
Cell type	adherent cells

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / case	20
Piece(s) / pallet	3360
Height of case	150 mm
Width of case	218 mm
Depth of case	288 mm
Volume of case	0.0094 cbm
Weight of case	1.52 kg
EAN inner box	4038917255479
EAN case	4038917255462

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT