



Product specification

DishRack, height: 370 mm, transparent, for 88 petri dishes with 92 mm Ø



Product information

Product number	93.1646
----------------	---------

Product descriptions

Name	DishRack, height: 370 mm, transparent, for 88 petri dishes with 92 mm Ø
Product number	93.1646
Product description	DishRack, height: 370 mm, material: PC, transparent, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

Material & colours

Product material	Polycarbonate (PC)
Color of product	transparent

Size

Height of product	370 mm
-------------------	--------

Product characteristics

Storage capacity	88
------------------	----

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Packaging

Piece(s) / pallet	108
Height of case	190 mm
Width of case	193 mm
Depth of case	383 mm
Volume of case	0.014 cbm
Weight of case	1.35 kg
EAN inner box	4038917087353
EAN case	4038917030458

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT