



Product specification

DishRack, height: 240 mm, transparent, for 52 petri dishes with 92 mm Ø



Product information

Product number	93.1647
----------------	---------

Product descriptions

Name	DishRack, height: 240 mm, transparent, for 52 petri dishes with 92 mm Ø
Product number	93.1647
Product description	DishRack, height: 240 mm, material: PC, transparent, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

Material & colours

Product material	Polycarbonate (PC)
Color of product	transparent

Size

Height of product	240 mm
-------------------	--------

Product characteristics

Storage capacity	52
------------------	----

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Batched

no

Packaging

Piece(s) / pallet	162
Height of case	190 mm
Width of case	193 mm
Depth of case	253 mm
Volume of case	0.0093 cbm
Weight of case	0.8 kg
EAN inner box	4038917087360
EAN case	4038917030502

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT