



Product specification

lumox® multiwell, Cell culture plate, with foil base, 24 well, 20 piece(s)



Product information

Product number	94.6110.024
----------------	-------------

Product descriptions

Name	lumox® multiwell, Cell culture plate, with foil base, 24 well, 20 piece(s)
Product number	94.6110.024
Product description	lumox® multiwell, Cell culture plate, with foil base, 24 well, frame: black, frame: PS, for adherent cells, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Product material	lumox
Frame material	Polystyrene (PS)
Color of frame	black

Size

Width of product	85.48 mm
Length of product	127.76 mm
Height of product	23.4 mm
Foil thickness	50 µm

Product characteristics

Growth area	1.9 cm ²
Number of wells	24
Design	with foil base
Base shape	Flat base
Vessel type	Plate
Closure type	Slip-on lid
Cell type	adherent cells

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	20
Piece(s) / case	20
Piece(s) / pallet	1920
Height of case	150 mm
Width of case	303 mm
Depth of case	333 mm
Volume of case	0.0151 cbm
Weight of case	1.67 kg
EAN inner box	4038917255448
EAN case	4038917255431

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT