



Product specification

Rack, PC, format: 6 x 3, suitable for tubes, S-Monovette®

13 mm Ø



Product information

Product number	93.827.013
----------------	------------

Product descriptions

Name	Rack, PC, format: 6 x 3, suitable for tubes, S-Monovette® 13 mm Ø
Product number	93.827.013
Product description	Rack, material: PC, transparent, format: 6 x 3, (LxWxH): 137 x 70 x 40 mm, for 18 collection tubes, suitable for tubes, S-Monovette® 13 mm Ø, 1 piece(s)/case

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	70 mm
Length of product	137 mm
Height of product	40 mm

Product characteristics

Storage capacity	18
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	36
Piece(s) / pallet	2376
Height of inner box	138 mm
Width of inner box	42 mm
Depth of inner box	72 mm
Height of case	155 mm
Width of case	268 mm
Depth of case	438 mm
Volume of case	0.0182 cbm
Weight of case	4.73 kg
EAN inner box	4038917121675
EAN case	4038917017923

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT