



Product specification

Rack, PC, format: 12 x 4, suitable for tubes, S-Monovette®

11 mm Ø



Product information

Product number	93.838
----------------	--------

Product descriptions

Name	Rack, PC, format: 12 x 4, suitable for tubes, S-Monovette® 11 mm Ø
------	--------------------------------------------------------------------

Product number	93.838
----------------	--------

Product description	Rack, material: PC, transparent, format: 12 x 4, (LxWxH): 257 x 90 x 40 mm, for 48 collection tubes, suitable for tubes, S-Monovette® 11 mm Ø, 1 piece(s)/case
---------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------

Material & colours

Product material	Polycarbonate (PC)
------------------	--------------------

Color of product	transparent
------------------	-------------

Size

Width of product	90 mm
------------------	-------

Length of product	257 mm
-------------------	--------

Height of product	40 mm
-------------------	-------

Product characteristics

Storage capacity	48
------------------	----

Autoclavable (max. °C)

121 °C

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	40
Piece(s) / pallet	1200
Height of inner box	259 mm
Width of inner box	42 mm
Depth of inner box	92 mm
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	10.29 kg
EAN inner box	4038917084352
EAN case	4038917003216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT