




Product specification

Rack, PC, format: 10 x 2, suitable for tubes up to 26 mm Ø



Product information

Product number	93.841.100
----------------	------------

Product descriptions

Name	Rack, PC, format: 10 x 2, suitable for tubes up to 26 mm Ø
Product number	93.841.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 327 x 72 x 60 mm, for 20 collection tubes, suitable for tubes up to 26 mm Ø, 1 piece(s)/case

Material & colours

Product material	Polycarbonate (PC)
Color of product	transparent

Size

Width of product	72 mm
Length of product	327 mm
Height of product	60 mm

Product characteristics

Storage capacity	20
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	40
Piece(s) / pallet	800
Height of case	360 mm
Width of case	387 mm
Depth of case	584 mm
Volume of case	0.0814 cbm
Weight of case	11.69 kg
EAN inner box	4038917121705
EAN case	4038917362368

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT