



Product specification

**Rack, PC, format: 10 x 2, suitable for tubes, S-Monovette®
13 mm and 11 mm Ø**



Product information

Product number	93.846.100
----------------	------------

Product descriptions

Name	Rack, PC, format: 10 x 2, suitable for tubes, S-Monovette® 13 mm and 11 mm Ø
Product number	93.846.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 257 x 62 x 40 mm, for 20 collection tubes, suitable for tubes, S-Monovette® 13 mm and 11 mm Ø, 1 piece(s)/case

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	62 mm
Length of product	257 mm
Height of product	40 mm

Product characteristics

Storage capacity	20
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	60
Piece(s) / pallet	1800
Height of inner box	260 mm
Width of inner box	43 mm
Depth of inner box	65 mm
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	11.99 kg
EAN inner box	4038917121729
EAN case	4038917017978

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT