



## Product specification

**Rack, PC, format: 6 x 3, suitable for tubes, S-Monovette®  
15 mm and 13 mm Ø**



## Product information

Product number	93.850
----------------	--------

## Product descriptions

Name	Rack, PC, format: 6 x 3, suitable for tubes, S-Monovette® 15 mm and 13 mm Ø
------	---

Product number	93.850
----------------	--------

Product description	Rack, material: PC, transparent, format: 6 x 3, (LxWxH): 137 x 70 x 40 mm, for 18 collection tubes, suitable for tubes, S-Monovette® 15 mm and 13 mm Ø, 1 piece(s)/case
---------------------	---

## Material & colours

Product material	Polycarbonate (PC)
------------------	--------------------

Color of product	transparent
------------------	-------------

## Size

Width of product	70 mm
------------------	-------

Length of product	137 mm
-------------------	--------

Height of product	40 mm
-------------------	-------

## Product characteristics

Storage capacity	18
------------------	----

Autoclavable (max. °C)

121 °C

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / case	36
Piece(s) / pallet	2376
Height of inner box	138 mm
Width of inner box	42 mm
Depth of inner box	72 mm
Height of case	155 mm
Width of case	268 mm
Depth of case	438 mm
Volume of case	0.0182 cbm
Weight of case	4.55 kg
EAN inner box	4038917084369
EAN case	4038917003223

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT