



## Product specification

**Block Rack D17, Ø opening: 17 mm, 5 x 10, yellow, with handle**



## Product information

Product number	93.851.172
----------------	------------

## Product descriptions

Name	Block Rack D17, Ø opening: 17 mm, 5 x 10, yellow, with handle
Product number	93.851.172
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, yellow, with handle, material: PP

## Material & colours

Color of product	yellow
------------------	--------

## Size

Width of product	109 mm
Length of product	246 mm
Height of product	45 mm
Diameter of opening	17 mm

## Product characteristics

Number of wells	50
Storage capacity	50
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / case	12
Piece(s) / pallet	840
Height of case	165 mm
Width of case	258 mm
Depth of case	453 mm
Volume of case	0.0193 cbm
Weight of case	2.79 kg
EAN inner box	4038917121842
EAN case	4038917074988

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT