



## Product specification

### Block Rack D17, Ø opening: 17 mm, 5 x 10, yellow



## Product information

Product number	93.852.172
----------------	------------

## Product descriptions

Name	Block Rack D17, Ø opening: 17 mm, 5 x 10, yellow
Product number	93.852.172
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, yellow, material: PP

## Material & colours

Product material	Polypropylene (PP)
Color of product	yellow

## Size

Width of product	109 mm
Length of product	209 mm
Height of product	45 mm
Diameter of opening	17 mm

## Product characteristics

Number of wells	50
Label/Print	With print
Storage capacity	50
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / case	20
Piece(s) / pallet	760
Height of case	120 mm
Width of case	433 mm
Depth of case	568 mm
Volume of case	0.0295 cbm
Weight of case	4.62 kg
EAN inner box	4038917122078
EAN case	4038917018135

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT