



## Product specification

### Block Rack D13, Ø opening: 13 mm, 5 x 10, blue



## Product information

Product number	93.853.134
----------------	------------

## Product descriptions

Name	Block Rack D13, Ø opening: 13 mm, 5 x 10, blue
Product number	93.853.134
Product description	Block Rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, blue, material: recycled PP

## Material & colours

Color of product	blue
------------------	------

## Size

Width of product	96 mm
Length of product	180 mm
Height of product	45 mm
Diameter of opening	13 mm

## Product characteristics

Number of wells	50
Label/Print	With print
Storage capacity	50
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / case	16
Piece(s) / pallet	1440
Height of case	120 mm
Width of case	368 mm
Depth of case	398 mm
Volume of case	0.0176 cbm
Weight of case	3.16 kg
EAN inner box	4038917122559
EAN case	4038917074056

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT