



Product specification

Block Rack D13, Ø opening: 13 mm, 5 x 10, orange



Product information

| | |
|----------------|------------|
| Product number | 93.853.135 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | Block Rack D13, Ø opening: 13 mm, 5 x 10, orange |
| Product number | 93.853.135 |
| Product description | Block Rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, orange, material: recycled PP |

Material & colours

| | |
|------------------|--------|
| Color of product | orange |
|------------------|--------|

Size

| | |
|---------------------|--------|
| Width of product | 96 mm |
| Length of product | 180 mm |
| Height of product | 45 mm |
| Diameter of opening | 13 mm |

Product characteristics

| | |
|------------------------|------------|
| Number of wells | 50 |
| Label/Print | With print |
| Storage capacity | 50 |
| Autoclavable (max. °C) | 121 °C |

Purity & certification

| | |
|------------------|---------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Batched | no |

Packaging

| | |
|-------------------|---------------|
| Piece(s) / case | 16 |
| Piece(s) / pallet | 1440 |
| Height of case | 120 mm |
| Width of case | 368 mm |
| Depth of case | 398 mm |
| Volume of case | 0.0176 cbm |
| Weight of case | 3.16 kg |
| EAN inner box | 4038917122566 |
| EAN case | 4038917074063 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT