


## Product specification

**Double block rack D13, Ø opening: 13 mm, 10 x 10, grey, connected lengthwise**



### Product information

Order number	93.853.200
--------------	------------

### Product descriptions

Name	Double block rack D13, Ø opening: 13 mm, 10 x 10, grey, connected lengthwise
------	--

Product number	93.853.200
----------------	------------

Product description	Double block rack D13, for 100 collection tubes, Ø opening: 13 mm, format: 10 x 10, grey, connected lengthwise, material: PP
---------------------	--

### Product characteristics

Number of wells	100
-----------------	-----

Storage capacity	100
------------------	-----

Autoclavable (max. °C)	121 °C
------------------------	--------

### Size

Width of product	192 mm
------------------	--------

Length of product	180 mm
-------------------	--------

Height of product	45 mm
-------------------	-------

Diameter of opening	13 mm
---------------------	-------

### Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Color of product	grey
------------------	------

Accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	no

## Packaging

Piece(s) / pallet	720
Height of case	120 mm
Width of case	368 mm
Depth of case	398 mm
Volume of case	0.0176 cbm
Weight of case	3.27 kg
EAN inner box	4038917122573
EAN case	4038917074070

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

