



## Product specification

### Block Rack D12, Ø opening: 12 mm, 5 x 10, yellow



## Product information

Product number	93.854.122
----------------	------------

## Product descriptions

Name	Block Rack D12, Ø opening: 12 mm, 5 x 10, yellow
Product number	93.854.122
Product description	Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, yellow, material: PP

## Material & colours

Color of product	yellow
------------------	--------

## Size

Width of product	80 mm
Length of product	152 mm
Height of product	45 mm
Diameter of opening	12 mm

## Product characteristics

Number of wells	50
Label/Print	With print
Storage capacity	50
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / case	20
Piece(s) / pallet	1560
Height of case	173 mm
Width of case	355 mm
Depth of case	390 mm
Volume of case	0.024 cbm
Weight of case	2.8 kg
EAN inner box	4038917122658
EAN case	4038917018494

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT