




Product specification

Cryobox, 75 x 75 x 52 mm, format: 5 x 5, for 25 collection tubes



Product information

Product number	93.872.425
----------------	------------

Product descriptions

Name	Cryobox, 75 x 75 x 52 mm, format: 5 x 5, for 25 Collection tubes
------	--

Product number	93.872.425
----------------	------------

Product description	Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, yellow, slip-on lid with ventilation function, cap: transparent, (LxWxH): 75 x 75 x 52 mm, format: 5 x 5, for 25 collection tubes, for CryoPure tubes 1.2–2.0 ml internal and external thread, 5 piece(s)/bag
---------------------	---

Material & colours

Product material	Polycarbonate (PC)
------------------	--------------------

Color of product	yellow
------------------	--------

Color of cap	transparent
--------------	-------------

Size

Width of product	75 mm
------------------	-------

Length of product	75 mm
-------------------	-------

Height of product	52 mm
-------------------	-------

Product characteristics

Design	with numerical coding at each aperture
Graduation	no
Storage capacity	25
Cap	slip-on lid with ventilation function
Closure type	Slip-on lid

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / case	20
Piece(s) / pallet	2640
Height of case	502 mm
Width of case	143 mm
Depth of case	156 mm
Volume of case	0.0112 cbm
Weight of case	1.41 kg
EAN inner box	4038917122818
EAN case	4038917420174

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT