


Product specification

Cryobox, 132 x 132 x 53 mm, format: 9 x 9, for 81 collection tubes



Product information

Order number 93.873.481

Product descriptions

Name Cryobox, 132 x 132 x 53 mm, format: 9 x 9, for 81 Collection tubes

Product number 93.873.481

Product description Cryobox, with numerical coding at each aperture, material: PC, yellow, slip-on lid with ventilation option, cap: transparent, (LxWxH): 132 x 132 x 53 mm, format: 9 x 9, for 81 collection tubes, for CryoPure tubes 1.2–2.0 ml internal and external thread, 5

Product characteristics

Design with numerical coding at each aperture

Graduation no

Storage capacity 81

Cap slip-on lid with ventilation option

Type of cap Slip-on lid

Suitable for CryoPure tubes 1.2–2.0 ml internal and external thread

Size

Width of product 132 mm

Length of product 132 mm

Height of product 53 mm

Material & colours

Product material	Polycarbonate (PC)
Color of product	yellow
Color of cap	transparent

Purity & certification

Accordance with MDD or IVDD	no medical device IVD
Lot or serial number	yes

Packaging

Piece(s) / case	20
Piece(s) / pallet	720
Volume of case	0.0197 cbm
Weight of case	3.91 kg
EAN inner box	4038917122863
EAN case	4038917415095

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

