




Product specification

Storage box, slip-on lid, PP, format: 9 x 9, for 81 collection tubes



Product information

Product number	93.876.281
----------------	------------

Product descriptions

Name	Storage box, Slip-on lid, PP, Format: 9 x 9, for 81 Collection tubes
Product number	93.876.281
Product description	Storage box, slip-on lid, material: PP, pink, format: 9 x 9, for 81 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm, 5 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	pink
Closure material	Polypropylene (PP)
Color of cap	transparent

Size

Width of product	130 mm
Length of product	130 mm
Height of product	45 mm

Product characteristics

Storage capacity	81
Closure type	Slip-on lid
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	20
Piece(s) / case	20
Piece(s) / pallet	1680
Height of case	470 mm
Width of case	264 mm
Depth of case	132 mm
Volume of case	0.0163 cbm
Weight of case	3.17 kg
EAN inner box	4038917161053
EAN case	4038917415224

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT