



Product specification

Storage box, hinged lid, PP, format: 10 x 10, for 100 collection tubes



Product information

| | |
|----------------|------------|
| Product number | 93.877.210 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | Storage box, Hinged lid, PP, Format: 10 x 10, for 100 Collection tubes |
| Product number | 93.877.210 |
| Product description | Storage box, hinged lid, material: PP, pink, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm, 5 piece(s)/bag |

Material & colours

| | |
|------------------|--------------------|
| Color of product | pink |
| Closure material | Polypropylene (PP) |
| Color of cap | pink |

Size

| | |
|-------------------|--------|
| Width of product | 140 mm |
| Length of product | 150 mm |
| Height of product | 55 mm |

| | |
|------------------------|------------|
| Storage capacity | 100 |
| Closure type | Hinged lid |
| Autoclavable (max. °C) | 121 °C |

Purity & certification

| | |
|------------------|---------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Batched | yes |

Packaging

| | |
|---------------------|---------------|
| Piece(s) / case | 20 |
| Piece(s) / pallet | 720 |
| Height of inner box | 140 mm |
| Width of inner box | 140 mm |
| Depth of inner box | 277 mm |
| Height of case | 305 mm |
| Width of case | 279 mm |
| Depth of case | 305 mm |
| Volume of case | 0.026 cbm |
| Weight of case | 4.1 kg |
| EAN inner box | 4038917161176 |
| EAN case | 4038917161169 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT