



## Product specification

**Storage box, hinged lid, PP, format: 10 x 10, for 100 collection tubes**



### Product information

Product number	93.877.410
----------------	------------

### Product descriptions

Name	Storage box, Hinged lid, PP, Format: 10 x 10, for 100 Collection tubes
Product number	93.877.410
Product description	Storage box, hinged lid, material: PP, yellow, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm, 5 piece(s)/bag

### Material & colours

Color of product	yellow
Closure material	Polypropylene (PP)
Color of cap	yellow

### Size

Width of product	140 mm
Length of product	150 mm
Height of product	55 mm

Storage capacity	100
Closure type	Hinged lid
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / case	20
Piece(s) / pallet	720
Height of inner box	140 mm
Width of inner box	140 mm
Depth of inner box	277 mm
Height of case	165 mm
Width of case	290 mm
Depth of case	610 mm
Volume of case	0.0292 cbm
Weight of case	4.1 kg
EAN inner box	4038917161213
EAN case	4038917161206

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT