



Product specification

Storage box, hinged lid, PP, format: 10 x 10, for 100 collection tubes



Product information

Product number	93.877.999
----------------	------------

Product descriptions

Name	Storage box, Hinged lid, PP, Format: 10 x 10, for 100 Collection tubes
Product number	93.877.999
Product description	Storage box, hinged lid, material: PP, colour mix, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm, 5 piece(s)/bag

Material & colours

Color of product	Colour mix
Closure material	Polypropylene (PP)
Color of cap	Colour mix

Size

Width of product	140 mm
Length of product	150 mm
Height of product	55 mm

Storage capacity	100
Closure type	Hinged lid
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / case	20
Piece(s) / pallet	720
Volume of case	0.026 cbm
Weight of case	4.09 kg
EAN inner box	4038917240376
EAN case	4038917240369

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT