



Product specification

Rack, PC, format: 10 x 2, suitable for tubes 100 x 21.5 mm



Product information

Product number	93.893.100
----------------	------------

Product descriptions

Name	Rack, PC, format: 10 x 2, suitable for tubes 100 x 21.5 mm
Product number	93.893.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 327 x 72 x 60 mm, for 20 collection tubes, suitable for tubes 100 x 21.5 mm, 1 piece(s)/case

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	72 mm
Length of product	327 mm
Height of product	60 mm

Product characteristics

Storage capacity	20
Autoclavable (max. °C)	121 °C

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	40
Piece(s) / pallet	800
Height of case	350 mm
Width of case	387 mm
Depth of case	584 mm
Volume of case	0.0791 cbm
Weight of case	11.57 kg
EAN inner box	4038917123136
EAN case	4038917362375

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT