



Product specification

CryoRack 40, PP, format: 10 x 4, suitable for CryoPure tubes



Product information

Product number	93.856.040
----------------	------------

Product descriptions

Name	CryoRack 40, PP, format: 10 x 4, suitable for CryoPure tubes
Product number	93.856.040
Product description	CryoRack 40, material: PP, format: 10 x 4, (LxWxH): 100 x 210 x 22 mm, for 40 collection tubes, suitable for CryoPure tubes, 1 piece(s)/bag

Size

Width of product	210 mm
Length of product	100 mm
Height of product	22 mm

Product characteristics

Graduation	no
Storage capacity	40

Purity & certification

Product category	not regulated
CE certified	no CE mark

Batched	no
---------	----

Packaging

Piece(s) / case	10
Piece(s) / pallet	1650
Height of inner box	25 mm
Width of inner box	120 mm
Depth of inner box	270 mm
Height of case	155 mm
Width of case	218 mm
Depth of case	228 mm
Volume of case	0.0077 cbm
Weight of case	0.9 kg
EAN inner box	4038917122788
EAN case	4038917077293

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT