




Product specification

lumox® multiwell, Cell culture plate, with foil base, 96 well, 4 piece(s)



Product information

Product number	94.6000.024
----------------	-------------

Product descriptions

Name	lumox® multiwell, Cell culture plate, with foil base, 96 well, 4 piece(s)
------	---

Product number	94.6000.024
----------------	-------------

Product description	lumox® multiwell, Cell culture plate, with foil base, 96 well, frame: black, frame: PS, for adherent cells, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister
---------------------	---

Material & colours

Frame material	Polystyrene (PS)
----------------	------------------

Color of frame	black
----------------	-------

Size

Width of product	85.48 mm
------------------	----------

Length of product	127.76 mm
-------------------	-----------

Height of product	17.2 mm
-------------------	---------

Foil thickness	50 µm
----------------	-------

Product characteristics

Growth area	0.34 cm ²
Number of wells	96
Design	with foil base
Base shape	Flat base
Vessel type	Plate
Closure type	Slip-on lid
Cell type	adherent cells

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	16
Piece(s) / pallet	1024
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	1.86 kg
EAN inner box	4038917254601
EAN case	4038917254595

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT