



Product specification

flexiPERM®, micro 12, re-usable silicon insert, 12 well, 1 piece(s)/bag



Product information

Product number	94.6011.416
----------------	-------------

Product descriptions

Name	flexiPERM®, micro 12, re-usable silicon insert, 12 well, 1 piece(s)/bag
Product number	94.6011.416
Product description	flexiPERM®, micro 12, re-usable silicon insert, for dividing the growth surface, 12 well, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Product material	silicone
------------------	----------

Size

Width of product	25 mm
Length of product	68 mm
Height of product	10 mm

Product characteristics

Growth area	0.3 cm ²
Number of wells	12

Purity & certification

Product category	no medical device IVD
Purity standard	non-pyrogenic/endotoxin-free, non-cytotoxic
Batched	yes

Packaging

Piece(s) / case	20
Piece(s) / pallet	1440
Height of case	225 mm
Width of case	253 mm
Depth of case	388 mm
Volume of case	0.0221 cbm
Weight of case	1.33 kg
EAN inner box	4038917267939
EAN case	4038917267922

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT