



Product specification

miniPERM®, HDC 50 bioreactor, 2-compartment bioreactor, for suspension cells



Product information

Product number	94.6077.121
----------------	-------------

Product descriptions

Name	miniPERM®, HDC 50 bioreactor, 2-compartment bioreactor, for suspension cells
Product number	94.6077.121
Product description	miniPERM®, HDC 50 bioreactor, 2-compartment bioreactor, for protein, virus and biomass production, for suspension cells, volume of product module: 50 ml, not autoclavable, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Size

Diameter	102 mm
Height of product	157 mm
Sample volume	400 ml
Production module volume	50 ml

Product characteristics

Design	for suspension cells
Closure type	Screw cap
Cell type	suspension cells

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / case	12
Piece(s) / pallet	240
Height of case	310 mm
Width of case	364 mm
Depth of case	498 mm
Volume of case	0.0562 cbm
Weight of case	5.4 kg
EAN inner box	4038917259132
EAN case	4038917259125

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT