



Product specification

lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, suspension cells



Product information

Product number	94.6077.305
----------------	-------------

Product descriptions

Name	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, suspension cells
------	---

Product number	94.6077.305
----------------	-------------

Product description	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, for suspension cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag
---------------------	---

Material & colours

Frame material	Polystyrene (PS)
----------------	------------------

Size

Diameter	50 mm
Height of product	12 mm
Foil thickness	25 µm

Product characteristics

Growth area	20.35 cm²
Design	with foil base

Base shape	Flat base
Vessel type	Dish
Closure type	Slip-on lid
Cell type	suspension cells

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	200
Piece(s) / pallet	7200
Height of case	260 mm
Width of case	337 mm
Depth of case	365 mm
Volume of case	0.032 cbm
Weight of case	3.52 kg
EAN inner box	4038917256070
EAN case	4038917256063

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT