



Product specification

Storage box, slip-on lid, PP, format: 4 x 4, for 16 collection tubes



Product information

Product number	93.881
----------------	--------

Product descriptions

Name	Storage box, Slip-on lid, PP, Format: 4 x 4, for 16 Collection tubes
------	--

Product number	93.881
----------------	--------

Product description	Storage box, slip-on lid, material: PP, blue, format: 4 x 4, for 16 collection tubes, suitable for 50 ml tubes and tubes up to Ø 30 mm, 2 piece(s)/bag
---------------------	--

Material & colours

Color of product	blue
------------------	------

Closure material	Polypropylene (PP)
------------------	--------------------

Color of cap	transparent
--------------	-------------

Size

Width of product	144 mm
------------------	--------

Length of product	144 mm
-------------------	--------

Height of product	127 mm
-------------------	--------

Product characteristics

Storage capacity	16
------------------	----

Closure type	Slip-on lid
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	12
Piece(s) / pallet	420
Height of inner box	145 mm
Width of inner box	145 mm
Depth of inner box	255 mm
Height of case	300 mm
Width of case	251 mm
Depth of case	452 mm
Volume of case	0.034 cbm
Weight of case	3.2 kg
EAN inner box	4038917355223
EAN case	4038917355216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT