




Product specification

**x-well cell culture chamber, 1 well, on lumox® slide,
removable frame**



Product information

Product number	94.6150.101
----------------	-------------

Product descriptions

Name	x-well cell culture chamber, 1 well, on lumox® slide, removable frame
Product number	94.6150.101
Product description	x-well cell culture chamber, 1 well, on lumox® slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

Material & colours

Product material	lumox
Frame material	Polystyrene (PS)

Size

Width of product	24 mm
Length of product	76 mm
Foil thickness	50 µm

Product characteristics

Growth area	9 cm ²
Removable frame	yes

Design	on lumox slide
Closure type	Slip-on lid
Cell type	adherent cells

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	96
Piece(s) / pallet	10368
Height of inner box	258 mm
Width of inner box	19 mm
Depth of inner box	127 mm
Height of case	185 mm
Width of case	270 mm
Depth of case	275 mm
Volume of case	0.0137 cbm
Weight of case	5.3 kg
EAN inner box	4038917255561
EAN case	4038917255554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT