




Product specification

**x-well cell culture chamber, 2 wells, on lumox® slide,
removable frame**



Product information

Product number	94.6150.201
----------------	-------------

Product descriptions

Name	x-well cell culture chamber, 2 wells, on lumox® slide, removable frame
Product number	94.6150.201
Product description	x-well cell culture chamber, 2 wells, on lumox® slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

Material & colours

Frame material	Polystyrene (PS)
----------------	------------------

Size

Width of product	24 mm
Length of product	76 mm
Foil thickness	50 µm

Product characteristics

Growth area	4.4 cm ²
Removable frame	yes
Design	on lumox slide

Closure type

Slip-on lid

Cell type

adherent cells

Purity & certification

Product category

not regulated

CE certified

no CE mark

Purity standard

sterile, non-pyrogenic/endotoxin-free,
non-cytotoxic

Sterilisation

Gamma irradiation

Pyrogen-free

yes

Batched

yes

Packaging

Piece(s) / case

96

Piece(s) / pallet

10368

Height of inner box

258 mm

Width of inner box

19 mm

Depth of inner box

127 mm

Height of case

185 mm

Width of case

270 mm

Depth of case

275 mm

Volume of case

0.0137 cbm

Weight of case

5.34 kg

EAN inner box

4038917255592

EAN case

4038917255585

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT