




# Product specification

**Storage box, slip-on lid, PP, format: 6 x 6, for 36 collection tubes**



## Product information

Product number	93.880
----------------	--------

## Product descriptions

Name	Storage box, Slip-on lid, PP, Format: 6 x 6, for 36 Collection tubes
Product number	93.880
Product description	Storage box, slip-on lid, material: PP, green, format: 6 x 6, for 36 collection tubes, suitable for 15 ml tubes and tubes up to Ø 17 mm, 2 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	green
Closure material	Polypropylene (PP)
Color of cap	transparent

## Size

Width of product	144 mm
Length of product	144 mm
Height of product	127 mm

## Product characteristics

Storage capacity	36
Closure type	Slip-on lid
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / case	12
Piece(s) / pallet	420
Height of case	452 mm
Width of case	300 mm
Depth of case	251 mm
Volume of case	0.034 cbm
Weight of case	3.2 kg
EAN inner box	4038917355193
EAN case	4038917355186

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT