




Product specification

**x-well cell culture chamber, 8 wells, on lumox® slide,
removable frame**



Product information

Product number	94.6150.801
----------------	-------------

Product descriptions

Name	x-well cell culture chamber, 8 wells, on lumox® slide, removable frame
------	--

Product number	94.6150.801
----------------	-------------

Product description	x-well cell culture chamber, 8 wells, on lumox® slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister
---------------------	--

Material & colours

Frame material	Polystyrene (PS)
----------------	------------------

Size

Width of product	24 mm
------------------	-------

Length of product	76 mm
-------------------	-------

Product characteristics

Growth area	0.8 cm ²
-------------	---------------------

Removable frame	yes
-----------------	-----

Design	on lumox slide
--------	----------------

Closure type	Slip-on lid
--------------	-------------

Cell type	adherent cells
-----------	----------------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	96
Piece(s) / pallet	10368
Height of inner box	258 mm
Width of inner box	19 mm
Depth of inner box	127 mm
Height of case	185 mm
Width of case	270 mm
Depth of case	275 mm
Volume of case	0.0137 cbm
Weight of case	5.48 kg
EAN inner box	4038917256223
EAN case	4038917256216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT