




Product specification

**x-well cell culture chamber, flask, on glass slide,
removable frame**



Product information

Product number	94.6170.002
----------------	-------------

Product descriptions

Name	x-well cell culture chamber, Flask, on glass slide, removable frame
Product number	94.6170.002
Product description	x-well cell culture chamber, flask, on glass slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

Material & colours

Product material	Soda-lime glass
Frame material	Polystyrene (PS)

Size

Width of product	24 mm
Length of product	76 mm

Product characteristics

Growth area	9 cm ²
Removable frame	yes
Design	on glass slide

Closure type	Screw cap
Cell type	adherent cells

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	96
Piece(s) / pallet	8064
Height of inner box	259 mm
Width of inner box	28 mm
Depth of inner box	127 mm
Height of case	255 mm
Width of case	270 mm
Depth of case	275 mm
Volume of case	0.0189 cbm
Weight of case	2.45 kg
EAN inner box	4038917254366
EAN case	4038917254359

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT