



Product specification

Rack, PP, format: 12 x 8, suitable for micro tubes and tubes with Ø 4-30 mm



Product information

Product number	93.879
----------------	--------

Product descriptions

Name	Rack, PP, format: 12 x 8, suitable for micro tubes and tubes with Ø 4-30 mm
Product number	93.879
Product description	Rack, material: PP, multi-coloured, format: 12 x 8, (LxWxH): 108 x 108 x 257 mm, suitable for micro tubes and tubes with Ø 4-30 mm, 1 piece(s)/bag

Material & colours

Color of product	multi-coloured
------------------	----------------

Size

Width of product	108 mm
Length of product	108 mm
Height of product	257 mm

Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / pallet	384
Height of inner box	110 mm
Width of inner box	117 mm
Depth of inner box	258 mm
Height of case	125 mm
Width of case	270 mm
Depth of case	445 mm
Volume of case	0.015 cbm
Weight of case	3.6 kg
EAN inner box	4038917355162
EAN case	4038917355155

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT