




## Product specification

**Rack, PP, format: 12 x 8, suitable for micro tubes and tubes with Ø 4-30 mm**



## Product information

Product number	93.879
----------------	--------

## Product descriptions

Name	Rack, PP, format: 12 x 8, suitable for micro tubes and tubes with Ø 4-30 mm
Product number	93.879
Product description	Rack, material: PP, multi-coloured, format: 12 x 8, (LxWxH): 108 x 108 x 257 mm, suitable for micro tubes and tubes with Ø 4-30 mm, 1 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	multi-coloured

## Size

Width of product	108 mm
Length of product	108 mm
Height of product	257 mm

## Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Piece(s) / pallet	384
Height of case	241 mm
Width of case	254 mm
Depth of case	216 mm
Volume of case	0.0132 cbm
Weight of case	3.6 kg
EAN inner box	4038917355162
EAN case	4038917355155

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT