



Product specification

Storage box, hinged lid, PP, format: 10 x 10, for 100 collection tubes



Product information

Product number	93.877.610
----------------	------------

Product descriptions

Name	Storage box, Hinged lid, PP, Format: 10 x 10, for 100 Collection tubes
Product number	93.877.610
Product description	Storage box, hinged lid, material: PP, blue, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm, 5 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	blue
Closure material	Polypropylene (PP)
Color of cap	blue

Size

Width of product	140 mm
Length of product	150 mm
Height of product	55 mm

Product characteristics

Storage capacity	100
Closure type	Hinged lid
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	20
Piece(s) / case	20
Piece(s) / pallet	720
Height of case	305 mm
Width of case	279 mm
Depth of case	305 mm
Volume of case	0.026 cbm
Weight of case	4.1 kg
EAN inner box	4038917161251
EAN case	4038917161244

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT