



# Product specification

## Cold transport container, english



### Product information

Product number	95.1123
----------------	---------

### Product descriptions

Name	Cold transport container, english
Product number	95.1123
Product description	Cold transport container, english

### Material & colours

Product material	Acrylnitrile butadiene styrene (ABS)
Color of product	white

### Size

Width of product	71 mm
Length of product	140 mm
Height of product	170 mm

### Product characteristics

Label/Print	With print
Color of print/label	red
Suction insert	no

### Purity & certification

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD

Batched	yes
---------	-----

## Packaging

Piece(s) / pallet	648
Height of case	195 mm
Width of case	223 mm
Depth of case	308 mm
Volume of case	0.0134 cbm
Weight of case	2.3 kg
EAN inner box	4038917085403
EAN case	4038917004084

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT